



SPHEREx Test and Calibration Control and Archive Software: User Guide

September 22, 2021
v0.5

Sam Condon scondon@caltech.edu
Marco Viero mverno@caltech.edu

Introduction

The SPHEREx Test and Calibration Control and Archive Software provides a Python based toolset allowing a user to interface with the suite of SPHEREx spectral and focus calibration instrumentation. The toolset allows a user to manually control instruments in a measurement setup or specify fully automated measurement runs. This document provides a practical guide on how to use this toolset.

Contents

1	Spectral Calibration Control	1
1.1	GUI	1
1.1.1	Startup	1
1.1.2	Manual Tab	2
1.1.3	Automation Tab	2

1 Spectral Calibration Control

This section details the control of the SPHEREx spectral calibration instrumentation. As of version v0.5, control of spectral calibration instrumentation can only be achieved using a GUI application, though a command-line interface is planned to be developed in future versions. In v0.5, the following instruments are integrated into the control interface:

1. Oriel/Newport CS260 Monochromator
2. Edmund Optics High-Speed Filter Wheel (referred to as NDF wheel)
3. Stanford Research Systems Sr510 and Sr830 Lockin Amplifiers
4. Labjack U6 (for warm shutter control)
5. Thorlabs S410C, S122C, and ajk;fdlsdlfsa;kjdjklfsa; detectors

1.1 GUI

Control of the spectral calibration instrumentation can be accomplished via a GUI. This section details how to use the GUI to run manual and automated measurements with the spectral calibration instrument suite.



1.1.1 Startup

To launch the latest version of the spectral calibration GUI, navigate to the *SPHEREx-Calibration-Automation* directory within the *SPHEREx-Lab-Tools* repository and run the following command:

```
python spherex_calibration.py spectral
```

Upon successful startup, a window similar to what is seen in figure (1) will appear. Take note of the two main tabs provided by the interface, *Automation* and *Manual*. The purpose of each of these tabs is exactly as it sounds. *Automation* can be used to create and run automated measurements while *Manual* allows manual control and readout of the instruments present in the setup. The operation of each of these tabs is described in detail in the subsequent sections.




Figure 1: Spectral Calibration Control Main Window

1.1.2 Manual Tab

The manual tab allows the user to manually control the instruments in the spectral calibration setup. As of version v0.5, the instruments that are integrated into the control include

1.1.3 Automation Tab